

Lab #3 Hints

Question 1: The summary data sheet should include all of the data in your lab book.

The Cl₂ Stock concentration should be converted to mg/mL in order to calculate the Cl₂ dosage. In order to do this, you must also account for the fact that you diluted the stock solution before you analyzed it with the HACH kit!

3 mL stock into 300 mL bottle → diluted 100x

$$\therefore \text{Cl}_2 \text{ Stock (mg/mL)} = \text{HACH reading (mg/L)} \times 100 \text{ (dilution factor)} \div 1000 \text{ mL/L}$$

You then need to calculate an average using the stock measurements from all groups. This value should be approximately 0.2 mg/mL of Cl₂.

The “Cl₂ Added (mL)” is the amount of chlorine stock solution that each group added to their samples (i.e. the numbers that were on the board).

The “Cl₂ Residual (mg/L)” is the HACH kit reading for the sample multiplied by the dilution factor (i.e. many of you diluted the sample 5x or more in the HACH cells in order to get a reading that was in range).

$$\begin{aligned} \text{“Cl}_2 \text{ Dosage (mg/L)”} &= \text{Cl}_2 \text{ Added (mL)} \times \text{Cl}_2 \text{ Stock Conc (mg/mL)} \\ &\div \text{sample volume (i.e. 300 mL)} \times 1 \text{ L} / 1000 \text{ mL} \end{aligned}$$

You can then do the plot and visually determine the breakpoint dosage.

Question 2: Keep in mind that the ammonia-chlorine reaction should be predictable based on stoichiometry.

Questions 3 through 5: These are all straight-forward discussion questions worth the same number of marks each (refer to the marking scheme).